



BACKGROUND

Reagents used for molecular diagnostics are typically stored frozen or at least refrigerated to maintain viability. The proprietary NeuDry process developed by NeuMoDx offers a simple, low cost, and rapid process using a conventional conveyer oven to produce ready to use, ambient condition stable reagents capable of providing long open-shelf life for active reagents used in molecular diagnostic testing. The performance of reagents produced using the NeuDry process across wide variety of molecular diagnostic applications is presented here.

METHODS

Three types of reagents formulated and optimized for both long term storage and in-use shelf life were produced using the NeuDry process - 1) Nucleic acid isolation reagents including NeuMag™ affinity particles, lytic enzyme, and encapsulated DNA/RNA process controls, 2) PCR Master mix, and 3) RT-PCR Master mix. Performance of the NeuDry reagents was initially evaluated against the liquid version of these reagents. Efficacy of the reagents was further evaluated by testing across a variety of clinical matrices encompassing both DNA and RNA targets. Long-term stability and in-use shelf life was determined via accelerated and real time stability studies on fully automated NeuMoDx Molecular System.

DRYING PROCESS

PCR OR RT-PCR FORMULATION NUCLEIC ACID EXTRACTION FORMULATION



DISPENSE





DRY, 15-20 MIN







READY TO USE MOLECULAR DIAGNOSTIC REAGENTS WITH LONG TERM AMBIENT AND **ON-BOARD STORAGE FOR FULLY AUTOMATED SYSTEMS**

Betty Wu, Elizabeth Craig, Crystal Nguyen, Hui Lin Lee, Lijie Gong, Will Marshall, Mark Olson, and S Brahmasandra, NeuMoDx Molecular, Inc, Ann Arbor, MI

RESULTS

NeuDry extraction reagents demonstrated >24 months shelf life and 63 days in-use stability in real time stability in a variety of matrices, such as plasma, serum, whole blood, urine and swab specimens. The NeuDry PCR/RT-PCR reagents, for both qualitative and quantitative assays, provided excellent amplification relative fluorescence increase (>80%). The PCR Master Mix formulated for the detection of GBS provided 25 months shelf life when stored at 18-28°C and up to 63 days in-use shelf life in real time stability studies. For both DNA and HCV, have demonstrated over 12 months stability and up to 63 days in-use shelf life with no loss of activity in real time and/or accelerated stability studies. The NeuDry process has also proven to be easily applicable to master mixes using different Taq polymerases without substantial re-formulation. The entire drying process takes less than 30 min resulting in over 10,000 tests produced in a single 8-hour shift with two operators and one conveyor oven. In addition, the continuous workflow enables significantly easier automation (& reduced cost) compared to batch processes such as lyophilization process.

NEUMODX GBS TEST STRIP: REAL TIME **STABILITY STUDY, 25 MONTHS**



• NeuMoDx Group B Strep (GBS) Assay is a qualitative in vitro diagnostic test designed to detect Group B Streptococcus (GBS) DNA.

• Three lots of GBS Test Strips were stored at ambient temperature (~ 23°C) and 28°C at an average 40-45% humidity for 25 months.

 Stability of the NeuMoDx GBS Test Strips by testing aged strips on a NeuMoDx System using GBS spiked into Lim Broth at concentration of 1000 CFU/mL. • 27 replicates of GBS positive samples and 3 negatives tested for each lot of

NeuMoDx GBS Test Strip per storage condition. • Aged group performed equally well as unaged Control group based on positive detection rate of 100% demonstrated for both storage conditions.





Pre-Aging	Aı	mbient (N=6) 2)	28°C (N=62)				
In-Use, Days	Avg	SD	% Pos	Avg Ct	SD	% Pos		
0	32.9	0.59	100%	32.8	0.55	100%		
15	32.6	0.13	100%	33.1	0.60	100%		
32	32.5	0.11	100%	32.6	0.36	100%		
43	33.6	0.13	100%	32.7	0.37	100%		
63	32.8	0.52	100%	33.2	0.53	100%		
ΔCt	-0.1			0.4				

- Two lots NeuMoDx GBS Test Strips were pre-aged as part of real time stability study for 12 months.
- The test strips were removed from the packaging, and stored at "in-use" conditions - average temperature at 28°C & average of 40-45% humidity - for up to 63 days
- Activity of the strips accessed by testing of the aged strips on a NeuMoDx System using GBS spiked into Lim Broth at concentration of 1000 CFU/mL.
- 31 replicates of GBS positive samples and 3 negatives tested for each lot of GBS Test Strip per pre-storage condition.
- The test strips aged 63 days at "in-use" conditions performed equivalently to day O based on the positive detection rate and Δ Ct from time O.

63 days

presented in the figures on the left.)

for 25 months.

processed on a NeuMoDx Molecular System. target and RNA sample process control with low SD of Ct values.

CONCLUSION

NeuDry reagents demonstrated excellent nucleic acid isolation from several clinical matrices as well as efficient and robust PCR/RT-PCR amplification using the eluted DNA/RNA directly. In addition, the NeuDry reagents have demonstrated a long shelf life of 25 months with an in-use life of 63 days.

NEUMODX EXTRACTION PLATE STABILITY – REAL TIME, 25 MONTHS





- RNA Sample Process Controls.
- with average 40-45% humidity for 25 months
- previously presented in the figures on the left.)
- Molecular System using aged and unaged extraction plates.

NEUMODX EXTRACTION PLATE STABILITY – IN-USE, 63 DAYS AT 25 MONTH



 NeuMoDx Extraction Plate contains all the reagents required for DNA and RNA extraction, such as lytic enzymes and Neu magnetic beads, as well as DNA and

• Three lots of Extraction Plates stored at ambient temperature (~ 23°C) and 28°C

 Stability of the Extraction Plates was evaluated by testing of the aged Extraction Plates on the NeuMoDx System using NeuMoDx HCV Assay (GBS Assay Data

• Stability of the Extraction Plates for RNA assays was evaluated with a model RNA assay, NeuMoDx HCV Assay. HCV viral particles were spiked into human plasma at concentration of 60 IU/mL (27 replicates) and processed on a NeuMoDx

• Aged group performed equally well based on the 100% amplification of both HCV target and the RNA Sample Process Control with low SD of Ct values.

• Two lots NeuMoDx Extraction Plates were pre-aged as part of real time stability study

 The Extraction Plates were removed from the packaging, and stored at "in-use" conditions - average temperature at 31°C & average of 40-45% humidity - for up to

• Stability of the Extraction Plates was evaluated by testing of the aged Extraction Plates on the NeuMoDx System using NeuMoDx HCV Assay (GBS Assay Data previously

• HCV viral particles were spiked into human plasma at concentration of 60 IU/mL and

• Aged group performed equally well based on the 100% amplification of both HCV

NEUMODX HBV QUANT ASSAY TEST STRIP: REAL TIME STABILITY - 13 MONTHS





	NeuMoDx Test Strips Shelf Life, 13 Months										
	Storage Condition	Time O		3 Months		7 Months		13 Months			
		Avg	SD	Avg	SD	Avg	SD	Avg	SD	Slope	D
HBV Conc [Log10 IU/mL]	Ambient	5.48	0.01	5.49	0.02	5.53	0.07	5.53	0.04	0.0038	0.3
	28°C	5.48	0.01	5.46	0.05	5.57	0.03	5.58	0.01	0.0094	0.9
	Ambient	3.41	0.05	3.42	0.05	3.52	0.03	3.52	0.06	0.0100	1.0
	28°C	3.41	0.05	3.42	0.05	3.52	0.02	3.52	0.04	0.0098	0.9
	Ambient	1.22	0.02	1.28	0.05	1.28	0.04	1.31	0.07	0.0065	0.0
	28°C	1.22	0.02	1.28	0.03	1.28	0.04	1.31	0.13	0.0015	0.

• Three lots of NeuMoDx HBV Quant Assay Test Strips were stored at ambient (~ 23°C) and 28°C with average 40-45% humidity for 13 months.

• Activity of the test strips was evaluated by testing the aged test strips on a NeuMoDx System using human plasma spiked with HBV at concentrations of 5.5, 3.5 (10 replicates at each level) and 1.4 log10 IU/mL (20 replicates) per lot for each storage condition

• Linear regression analysis demonstrated minimal drift after 13 months - ≤ 1%, in quantification of HBV after storage for 13 months.

NEUMODX HBV QUANT ASSAY TEST STRIP: IN-USE STABILITY – 63 DAYS



Lot 1, 63 Days, (28°C Pre-aging Group) 1 3 5 7 9 11 13 15 17 19 21 23 25 27 29 31 33 35 37 39 41

Expected HBV Conc, [Log10 IU/mL]	DO		D15		D32		D43		D63		Δ fr DO	
	Avg	SD	mL]									
5.63	5.58	0.02	5.56	0.01	5.59	0.01	5.59	0.09	5.60	0.03	0.02	
3.63	3.51	0.01	3.47	0.01	3.54	0.00	3.52	0.05	3.55	0.05	0.03	
1.40	1.28	0.04	1.27	0.13	1.41	0.04	1.33	0.11	1.28	0.03	0.00	

• In-Use stability study of NeuMoDx HBV Test Strips. HBV Test Strips were first aged at ambient or 28°C as part of the Real Time Stability study. Two lots of HBV Test Strips were removed from their primary storage condition and stored at In-Use condition (31°C with average 40-45% humidity) up to 61 days.

• The activity of the strips was evaluated by testing the aged strips on the NeuMoDx System using human plasma spiked with HBV at concentrations of 5.5, 3.5 (10 replicates at each level) and 1.4 log10 IU/mL (20 replicates) per lot for each storage condition.

• There is minimum drift, ≤ 0.07%, in guantification of HBV up to 61 days at the In-Use condition for both pre-aging groups, ambient or 28°C. The data from the 28 °C pre-aging group are presented here.

ACKNOWLEDGMENTS

The authors gratefully acknowledge the help and support provided by all members of the NeuMoDx team. We would also like to thank our collaborators for their invaluable assistance in providing access to clinical specimens for testing and evaluation. NeuMoDx GBS Test and NeuMoDx Extraction Plates on NeuMoDx 288 is For In Vitro Diagnostics only. NeuMoDx HBV, CMV and HCV tests are For Research Use Only. Not for use in diagnostic procedures.



- Three lots of HCV Test Strips were stored at 40°C up to 149 days, equivalent to 17 months at 23°C. • The activity of the strips was evaluated by testing the aged strips on the NeuMoDx System using human
- plasma spiked with HCV at concentrations of 5, 3 and 1.88 log10 IU/mL (6 replicates per level) per lot for each storage condition.
- After storage up to 149 days at 40°C, the test strips quantified HCV targets with minimum bias to the expected concentrations. The difference between D149 and D0 is < 0.17 Log10 IU/mL and almost no drift, ≤ 0.07%, was observed when analyzed with linear regression.
- Three lots of test strips performed reproducibly with SD < 0.17 Log10 IU/mL across all the testing points and concentrations.
- Representative PCR amplification curves from one of the three lots of HCV Test Strips are shown above.

0.06%