

QUANTITATIVE DETECTION OF HCV USING THE NEUMODX MOLECULAR DIAGNOSTIC SYSTEM

Jessica Zhu*, Elizabeth Craig, Maureen Carey, Catherine Couture, Michelle Mastronardi, Betty Wu, NeuMoDx Molecular, Ann Arbor, MI

GENOTYPE SENSITIVITY

negative samples.

The LoD and LLoQ of the NeuMoDx HCV

negative plasma spiked with different HCV

genotypes at four different levels including

6.0

HCV genotype inclusivity. The NeuMoDx HCV

of HCV in negative plasma with LoD at or under

REFERENCE REFERENCE

POSITIVE AGREEMENT = 99.7% 95% CI (98.5% - 100%)

NEGATIVE AGREEMENT = 95.1% 95% CI (91.7% - 97.2%)

NeuMoDx HCV Test method correlation. The NeuMoDx HCV Test

demonstrates excellent correlation with the comparator HCV test

TESTS

POSITIVE

TABLE**

NMD_x TEST

NEGATIVE

results from a reference lab.

95%CI (-0.43,-0.14)

ope Coefficien

95%CI (0.99,1.05)

TESTS

NEGATIVE

285

the WHO 5th IS LoD.

Test accurately detected all six relevant genotypes

LoD (IU/mL) LLoQ (IU/mL)

Test was confirmed with pooled HCV

BACKGROUND Determining Hepatitis C Virus (HCV) RNA levels in plasma and/or serum is an important tool to characterize viral loads in infected patients to monitor disease progression, efficacy of antiviral therapies, as well as to detect drug resistant mutants and identify relapse upon discontinuation of an antiviral therapy. The NeuMoDx HCV Test is an in-vitro diagnostic assay incorporating a universal nucleic acid isolation chemistry enabling extraction of qPCR ready RNA from serum and plasma specimens, combined with a sensitive quantitative rt-PCR assay to deliver highly accurate results in a completely automated, "random access" manner on the NeuMoDx Molecular System. In addition, all reagents and disposables are room temperature stable and are intended to remain on-board the system to provide a seamless, on-demand testing workflow.


The NeuMoDx Molecular System automates and integrates the extraction, purification, quantification, quantification, quantification, and results interpretation of infectious disease nucleic acid targets using quantitative RT-PCR. The objective of this study was to test and report performance of the NeuMoDx HCV Test in key analytical performance metrics. Internal pre-analytical studies were performed to characterize the analytical sensitivity, linearity, precision, inclusivity, turnaround time, as well as characterizing quantitative correlation to a reference test using split samples. The results of these studies are presented here.

RESULTS The NeuMoDx HCV Test showed a detection limit of 7.5 IU/mL (95% CI of (6.4, 9.2)) and lower limit of quantification of 7.7 IU/mL using the 5th International WHO HCV Standard, with an overall LLoQ of 0.9 Log10 IU/mL for plasma and serum. A master calibration curve traceable to the 5th WHO HCV Standard as well as external calibrators (based on secondary standards traceable to the 5th WHO HCV Standard) were developed to provide accurate quantitative results across multiple systems and reagent lots. The NeuMoDx HCV Test showed equivalent detection performance across all relevant HCV genotypes and a time to first results of ~80 min. No cross-reactivity or interference was observed against any of the pathogens or agents tested. The method correlation study performed with the NeuMoDx HCV test and using split samples demonstrated excellent concordance with a reference test. An R² >0.95 and bias of HCV viral load output from the two tests of less than 0.5 Log IU/mL was obtained.

HCV Analytical Sensitivity (LoD & LLoQ)

PLASMA

The Limit of Detection of the NeuMoDx HCV Test was determined with pooled HCV negative plasma spiked with 5th WHO International Standard at seven different levels including negative samples. The limit of detection of HCV was determined to be 7.5 IU/mL based on Probit style analysis and the calculated LLoQ was determined to be 7.7 IU/mL.



Limit of Detection of the NeuMoDx HCV Test. Probit analysis from the data in the above table was used to determine the LoD of the HCV target to be 7.5IU/mL with a 95% CI of (6.4,9.2).

	NeuMoDx HCV Test LLoQ						
Target Conc. (IU/mL)	Target Conc. (Log ₁₀ IU/mL)	N	Abs. Bias	Standard Deviation (SD)	Total Analytical Error (TAE)	LLoQ (IU/mL)	
30.00	1.48	108	0.07	0.32	0.71		
15.00	1.18	108	0.06	0.36	0.79	7.7 IU/mL	
10.00	1.00	108	0.07	0.35	0.77	or	
7.50	0.88	108	0.13	0.44	1.02	0.89 Log ₁₀	
3.75	0.57	108	0.51	0.43	1.38	IU/mL	
1.88	0.27	108	0.83	0.36	1.55		

Lower limit of quantitation (LLoQ) of NeuMoDx HCV Test. The lowest target level detected at a rate > 95% AND with TAE (bias+2*SD) ≤ 1.0 was used to determine the LLoQ. The LLoQ of the NeuMoDx HCV Test was determined to be 7.7 IU/mL.

Method Correlation

The quantitative performance of the NeuMoDx HCV Test was

A total of 646 clinical plasma and serum specimens within

the linear range were used to generate the linear regression.

single-blinded study using clinical samples obtained from a

quantitative correlation with the reference test with Deming

Testing was performed internally at NeuMoDx through a

assessed against 4 CE/FDA approved comparator tests by


testing clinical specimens from HCV infected patients.

The NeuMoDx HCV test demonstrated excellent

and Passing-Bablok linear regression tests.

The Limit of Detection of the NeuMoDx HCV Test was determined with pooled HCV negative plasma spiked with AcroMetrix High Control HCV (traceable to the WHO 5th IS) at seven different levels including negative samples. The limit of detection of HCV was determined to be 8.0 IU/mL based on Probit style analysis and the calculated LLoQ was determined to be 8.4 IU/mL.

Ne	euMoDx HCV	Test	Limit of	Detection	on
Target Conc. (IU/mL)	Target Conc. (Log ₁₀ IU/mL)	N	# Positive	% Positive	LoD (Probit) (IU/mL)
30	1.48	108	108	100%	
15	1.18	108	107	99%	
10	1.00	108	102	94%	8.0 IU/mL
7.5	0.88	108	105	97%	95% CI (6.6-10.4
3.75	0.57	108	86	80%	`IU/mL)
1.875	0.27	108	63	58%	
NEG	N/A	107	1	0.93%	



Limit of Detection of the NeuMoDx HCV Test. Probit analysis from the data in the above table was used to determine the LoD of the HCV target to be 8.0IU/mL with a 95% CI of (6.6,10.4).

NeuMoDx HCV Test LLoQ						
Target Conc. (IU/mL)	Target Conc. (Log ₁₀ IU/mL)	N	Abs. Bias	Standard Deviation (SD)	Total Analytical Error (TAE)	LLoQ (IU/mL)
30.00	1.48	108	0.08	0.30	0.69	
15.00	1.18	108	0.06	0.32	0.70	8.4 IU/mL
10.00	1.00	108	0.14	0.36	0.85	or
7.50	0.88	108	0.25	0.25	1.09	0.9 Log ₁₀
3.75	0.57	108	0.59	0.58	1.76	IU/mL
1.88	0.27	108	0.84	0.69	2.22	

Lower limit of quantitation (LLoQ) of NeuMoDx HCV Test. The lowest target level detected at a rate > 95% AND with TAE (bias+2*SD) ≤ 1.0 was used to determine the LLoQ. The LLoQ of the NeuMoDx HCV Test was determined to be 8.4 IU/mL.

PASSING-BABLOCK:

Regression of NMDx Concentration by

Reference Test Concentration

Reference Test Concentration



NeuMoDx Molecular System Streamlined Testing

- specimens to providing real-time PCR results in
- True Random Access: Ability to mix specimen
- hour shift for the N288, ~100 RNA tests in an 8 hour shift for the N96
- Continuous Loading of Specimens: Specimens and Reagents can be loaded/unloaded at any time Large Walk-Away Window: Up to 288 samples for the N288, 96 samples for the N96 Seamless On Demand Operation: Automated inventory management of consumables and reagents
- Long In-Use Shelf Life: On-board room temperature stable reagents Real-time PCR: Five-color fluorescence det offers real-time PCR multiplexing ability






HCV Linearity

PLASMA

The linearity of the NeuMoDx HCV Test was determined by diluting either Acrometrix HCV RNA or Asuragen HCV Armored RNA® in pooled HCV negative plasma to create a panel spanning 10 logs of HCV concentration ranging from 8.0 Log10 IU/mL to 1.5 Log10 IU/mL. Data from this study showed that the NeuMoDx HCV Test demonstrated excellent linearity across the 10 logs.

	NeuMoDx HCV Test Linearity							
Target Conc. (IU/mL)	Target Conc. (Log ₁₀ IU/mL)	Mean Conc. (Log ₁₀ IU/mL)	Bias	Calculated Linear Fit				
1.00E+08	8.00	8.00	0.00	8.00				
1.00E+07	7.00	7.03	0.03	7.00				
1.00E+06	6.00	6.04	0.04	6.00				
3.16E+05	5.50	5.54	0.04	5.50				
1.00E+05	5.00	4.95	0.05	5.00				
3.15E+04	4.50	4.46	0.04	4.50				
1.00E+04	4.00	3.97	0.03	4.00				
3.16E+03	3.50	3.46	0.04	3.50				
3.16E+02	2.50	2.52	0.02	2.50				
3.16E+01	1.50	1.55	0.05	1.50				



- Linear range of the NeuMoDx HCV Test. The NeuMoDx HCV Test is linear over 10 Log10 units.
- Linear range of the NeuMoDx HCV Test. The NeuMoDx HCV Test shows excellent correlation over 10 Log10 units.

SERUM

Fast Time to First Results: ~75 min

The linearity of the NeuMoDx HCV Test was determined by diluting either Acrometrix HCV RNA or Asuragen HCV Armored RNA® in pooled HCV negative serum to create a panel spanning 8 logs of HCV concentration ranging from 8.0 Log10 IU/mL to 1.5 Log10 IU/mL. Data from this study showed that the NeuMoDx HCV Test demonstrated excellent linearity across the 8 logs.

NeuMoDx HCV Test Linearity						
Target Conc. (IU/mL)		Mean Conc. (Log ₁₀ IU/mL)	Bias	Calculated Linear Fit		
1.00E+08	8	8.54	0.54	8.48		
1.00E+07	7	7.54	0.54	7.43		
1.00E+06	6	6.46	0.46	6.38		
1.00E+05	5	5.41	0.41	5.33		
3.16E+04	4.5	4.77	0.27	4.80		
3.16E+03	3.5	3.73	0.22	3.75		
3.16E+02	2.5	2.74	0.24	2.70		
3.16E+01	1.5	1.71	0.21	1.65		

Linear range of the

NeuMoDx HCV Test.

Linear range of the

NeuMoDx HCV Test.


over 8 Log10 units.

The NeuMoDx HCV Test

shows excellent correlation

The NeuMoDx HCV Test is


linear over 10 Log10 units.



GENOTYPE SENSITIVITY

The linearity of the NeuMoDx HCV Test was determined by diluting each HCV genotype in pooled HCV negative plasma to create a panel spanning 6 logs of HCV concentration ranging from 6.5 Log10 IU/mL to 0.8 Log10 IU/mL (depending on concentration of each genotype). Data from this study showed that the NeuMoDx HCV Test demonstrated excellent linearity across the 6 logs.

HCV Genotype Panel	LoD (IU/mL)	LLoQ (IU/mL)
1	y = 1.054x + 0.1325	0.9790
2	y = 1.0792x - 0.0748	0.9854
3	y = 1.0423x - 0.0439	0.9814
4	y = 1.0158x + 0.0292	0.9729
5	y = 0.9873x + 0.1524	0.9939
6	y = 1.0393x + 0.0396	0.9973
		·



Cross-Reactivity & Interference

Adenovirus 2	Dengue virus 4	HIV-1	Influenza A virus	Zika virus	Staphylococcus epidermidis
Adenovirus 5	Epstein-Barr virus	HIV-2	Parvovirus B19	Candida albicans	
Cytomegalovirus	Hepatitis A Virus	Human papillomavirus 16	Rubella virus	Propionibacterium acnes	
Dengue virus 1	Hepatitis B Virus	Human papillomavirus 18	St. Louis encephalitis virus	Staphylococcus aureus	
Dengue virus 2	Human herpes virus 6B	Human T-cell lymphotropic virus-type 1	Yellow Fever	Neisseria gonorrhoeae	
Dengue virus 3	Human herpes virus 8	Human T-cell lymphotropic virus-type 2	West Nile	Chalymidia trachomatis	

Non-Target Organisms	Average conc. (Log10 IU/mL)	Bias* (Log10 IU/mL)
Pool 1: Cytomegalovirus, Epstein-Barr virus, Human Herpesvirus 6B, Human Herpesvirus 8, Rubella	1.71	0.31
Pool 3: Hepatitis A, Hepatitis B, Zika virus, Yellow Fever, Human Immunodeficiency virus 1	0.99	0.41
Pool 6: Candida albicans, Chlamydia trachomatis, Herpes Simplex virus 1, Herpes Simplex virus 2, Staphylococcus aureus, Staphylococcus epidermidis	1.44	0.04

Non-Target Organisms	Average conc. (Log10 IU/mL)	Bias* (Log10 IU/mL)
Adenovirus 2	1.75	0.35
Adenovirus 5	1.09	0.31
Dengue 4	1.61	0.21
Parvo B19	1.76	0.36
Dengue 2	1.55	0.15
Dengue 3	1.67	0.27
St Louis Encephalitis	1.64	0.25
Dengue 1	1.77	0.37
Human T-lymphotropic virus 1	1.59	0.19
Human T-lymphotropic virus 2	1.55	0.15
Influenza A	1.6	0.20
Human Immunodeficiency virus 2	1.54	0.15
Human papillomavirus 16	1.5	0.10
Human papillomavirus 18	1.53	0.13
Neisseria Gonorrhoea	1.64	0.24
Propionibacterium acnes	1.41	0.01
	1	

HCV Interfering Substances – Commensal Organisms. The same high titers of the same 31 non-target organisms (pooled or individual) tested. The NeuMoDx HCV Test had minimal deviation of quantitation from HCV control samples, with no significant interference.

West Nile

1.65 0.25

▲ HCV Analytical Specificity - Cross-Reactivity. The NeuMoDx HCV Test had no cross-reactivity with 31 non-target (or phylogenetically similar) organisms found in blood/plasma specimens – demonstrated 100% analytical specificity.

Average conc. Bias*
(Log10 IU/mL) (Log10 IU/mL)

Bilirubin	1.47	0.14
Protein (albumin)	1.47	0.14
Hemoglobin	1.61	0.28
Triglycerides	1.31	0.02
Disease State	Average conc. (Log10 IU/mL)	Bias* (Log10 IU/mL)
Systemic Lupus Erythematosus (SLE)	1.29	0.06
Antinuclear Antibody (ANA)	1.53	0.18
Rheumatoid Arthritis (RA)	1.39	0.04
Rheumatoid Factor (RF)	1.43	0.08
Non-Alcoholic Steatohepatitis (NASH)	1.32	0.03
HBV Antibodies	1.45	0.10
HCV Antibodies	1.43	0.08
Alcoholic cirrhosis (AC)	1.29	0.06
Non Townst Overniems	Average con	c. Bias*

Non-Target Organisms	Average conc. (Log10 IU/mL)	Bias* (Log10 IU/mL)
Pool 1: Sofosbuvir, Ledipasvir, Velpatasivr, Clarithromycin, Interferon alfa-2a	1.48	0.15
Pool 2: Paritaprevir, Ombitasvir, Ritonavir, Abacavir sulfate, Ribavirin	1.4	0.07
Pool 3: Grazoprevir, Elbasvir, Tenofovir disoproxil, Lamivudine, Valganciclovir	1.4	0.07
Pool 4: Efavirenz, Lopinavir, Azithromycin, Dolutegravir, Simeprevir	1.51	0.18
Pool 5: Emtricitabine, Raltegravir, Amoxicillin, Rilpivirine, Dasabuvir, Glecaprevir	1.4	0.07

HCV Interfering Substances – Endogenous & Exogenous Substances. Twelve (12) different endogenous/disease state specimens & 26 exogenous substances (drugs) tested. The NeuMoDx HCV Test had minimal deviation of quantitation from HCV control samples, with no significant interference.

Precision

The Within Lab Precision of the NeuMoDx HCV test was determined by testing a 7 member panel of HCV on multiple NeuMoDx Molecular Systems across multiple days.

The same concentrations of HCV were used across these

consumables and compared to each other for consistency. The precision Within-run and Across-runs was characterized and the standard deviation for both was determined to be ≤ 0.5 Log10 IU/mL.

Precision of the NeuMoDx HCV Test.

The NeuMoDx HCV Test demonstrated excellent within laboratory precision calculated from the quantitative data across target levels, reagent lots, and Systems with a maximum overall standard deviation ≤ 0.26 Log10 IU/mL.

	NeuN	oDx HC\	/ Test	Within L	ab Prec	ision	
Panel Member	Target Conc. (Log10 IU/mL)	Mean Conc. (Log10 IU/mL)	N	Within System SD	Within Day SD	Within Run SD	(Overall) Within Lab SD
1	6.0	5.95	216	0.17	0.13	0.10	0.17
2	5.0	4.87	216	0.20	0.14	0.12	0.20
3	3.0	2.89	217	0.19	0.17	0.17	0.19
4	4.4	4.45	216	0.12	0.10	0.08	0.13
5	3.4	3.45	216	0.12	0.12	0.11	0.13
6	2.4	2.41	220	0.17	0.15	0.15	0.17
7	1.4	1.40	217	0.26	0.25	0.25	0.24

CONCLUSIONS & ACKNOWLEDGMENTS

DEMING:

Regression of NMDx Concentration by

Reference Test Concentration

Reference Test Concentration

The performance data presented here clearly demonstrate that the NeuMoDx HCV Test is an extremely easy to use, rapid, automated test for the sensitive and accurate monitoring of HCV viral load. The authors gratefully acknowledge the help and support provided by all members of the NeuMoDx team. We would also like to thank our collaborators for their invaluable assistance in providing access to clinical specimens for testing and evaluation.

Deming Analysis

95%CI (-0.37,0.06)

95%CI (0.966,1.03)

reference laboratory.