

# ANALYTICAL PERFORMANCE OF THE NEUMODX CTNG TEST ON THE NEUMODX MOLECULAR SYSTEMS FOR DETECTION OF CHLAMYDIA TRACHOMATIS AND NEISSERIA GONORRHOEAE FROM URINE AND SWAB SPECIMENS B Wu<sup>1</sup>, J Zhu<sup>1</sup>, M Carey<sup>1</sup>\*, C Lounds<sup>1</sup>, A Narwold<sup>1</sup>, M Mastronardi<sup>1</sup>, D Tomlinson, S Brahmasandra<sup>1</sup> <sup>1</sup>NeuMoDx Molecular, Ann Arbor, MI United States <sup>2</sup>Qiagen, Manchester, United Kingdom

## INTRODUCTION

Laboratories around the world rely on automated nucleic acid testing systems (NAATs) for screening and diagnosis of Chlamydia trachomatis (CT) and Neisseria gonorrhoeae (NG) infections in symptomatic/asymptomatic patients from urine and swab matrices. The NeuMoDx CTNG Test is CEIVD certified for use with urine specimens on the high-throughput floor model NeuMoDx 288 Molecular System (N288). The goal of this study was to demonstrate performance of the NeuMoDx CTNG Test on the medium throughput, bench-top NeuMoDx 96 System as well as across additional specimen types such as vaginal and endocervical swabs. Results of verification of these studies in both urine and swab matrices on both NeuMoDx Platforms are presented here.

## RESULTS

The NeuMoDx CTNG Test demonstrated a Limit of Detection of 9 EB/mL for CT and 5 cells/mL for CT and 5 cells/mL for NG in urine, and 20 EB/mL for CT and 5 cells/mL for NG in swabs with excellent inclusivity across the 15 CT serovars and ~10 NG isolates. The sensitivity of each target (CT or NG) was maintained in the presence of high target concentration (>1E6 cells/mL or EB/mL) of the other target (NG or CT) demonstrating the robustness of the assay. 100% specificity was observed across ~113 phylogenetically similar or co-habiting pathogens; and the test demonstrated no adverse effect in the presence of ~150 endogenous/ exogenous interfering agents or pathogens. The NeuMoDx CTNG Test demonstrated excellent reproducibility across multiple reagent lots. Time to first results was only 60 minutes enabling delivery of rapid results for emergency samples and the specimens were stable for at least 24 hours on board the worktable of the instrument enabling walkaway operation.

### Analytical Sensitivity (LoD)

The Limit of Detection of the NeuMoDx CTNG Test in the urine matrix was determined with pooled negative clinical urine spiked with AcroMetrix<sup>TM</sup> CT or NG at six different levels including negative samples.

The Limit of Detection of the NeuMoDx CTNG Test in swab matrices was determined in a hit-rate style with pooled negative clinical vaginal or endocervical swabs spiked with AcroMetrixTM CT or NG at three different levels including negative samples The limit of detection in urine was determined to be 6 EB/mL for CT and 5 cells/mL for NG (after interference studies, as well). The limit of detection in swabs was determined to be 20 EB/mL for CT and 5 cells/mL for NG.







overall LoD of the CT target in swab was determined to be 20 EB/mL.

 
 Table 1. Limit of Detection of the NeuMoDx CTNG Test in the urine matrix for the CT target.
 The overall LoD of the CT target in urine was determined to be 6 EB/mL.

Log[CONC/mL]

\*N96 instrument LoD was deemed acceptable due to smaller sample size

NeuMoDx CTNG Test Limit of Detection – CT (N288 & N96)								
Target Conc. (EB/mL)	Ν	# Positive	% Positive	LoD (Probit)				
	Vaginal Swab Results							
30	48	48	100%					
20	48	48	100%	20 FR/ml				
0	48	0	0%	LD/IIIL				
Endocervical Swab Results								
30	48	48	100%					
20	48	48	100%	20 FR/ml				
0	48	0	0%	ED/ML				

Log[CONC/mL]

 
 Table 2.
 Limit of Detection of the NeuMoDx CTNG Test in
 the swab matrix for the CT target. The overall LoD of the CT target in swab was determined to be 20 EB/mL.

NeuMoDx CTNG Test Limit of Detection								
	– NG (N288 & N96)							
Target Conc. (cells/mL)	Ν	# Positive	% Positive	LoD (Detection)				
	Vaginal Swab Results							
9	48	48	100%	_				
5	48	47	98%	5 cells/ml				
0	48	0	0%	Cens/IIIE				
Endocervical Swab Results								
9	48	48	100%	_				
5	48	48	100%	5 cells/ml				
0	48	0	0%	CEII3/IIIE				

 
 Table 4.
 Limit of Detection of the NeuMoDx CTNG Test in
 the swab matrix for the NG target. The overall LoD of the NG target in swab was determined to be 5 cells/mL.

CT Serovar	NG Clinical Isolate (ATCC #)
Α	49981
В	31426
Ba	31407
С	27633
LGV I	9793
LGV II	43070
LGV III	51109
E	35542
F	35541
G	49498
н	49926
I	
•	<b>Table 5</b> Limit of Detection of the NeuMoDy

**Table 5.** Limit of Defection of the NeulVioDx CTNG Test in the swab matrix for the NG target. The overall LoD of the NG target in swab was determined to be 5 cells/mL.

Analytical studies were performed in urine and both swab matrices across

the NeuMoDx 288 and NeuMoDx 96 systems to characterize various analytical parameters such as – sensitivity, cross-reactivity, inclusivity, impact of interfering substances, cross-contamination, specimen stability, and turnaround time (TAT).

0.25 cells/mL (95% CI (0.25, 0.25) cells/mL)

0.3 cells/mL (95% CI (0.28, 0.28) cells/mL)

### Table 3. Limit of Detection of the NeuMoDx CTNG Test in the swab matrix for the CT target. The

## Analytical Specificity & Interference

### **Analytical Specificity** & Cross-Reactivity

A total of 113 culture isolates or DNA from organisms potentially cohabiting or phylogenetically similar to either CT or NG were evaluated for possible crossreactivity with the NeuMoDx CTNG Test. Most of the organisms were spiked into CTNG negative urine at ~1E6 CFU/mL, except in the case of a few organisms where high copies of DNA (10 ng/mL) were spiked into CTNG negative urine. The list of organisms tested is shown in Table 6. In addition, it was demonstrated that the detection of CT and NG targets, even at low levels of target (~3X LoD), was not affected by the presence of any of these organisms

### Interfering Substances

The performance of the NeuMoDx CTNG Test in the presence of endogenous/ exogenous substances was evaluated by testing for efficacy of detection of low (3x LoD) and moderate (10x LoD) level CT & NG targets. This was done in both negative clinical urine and swab matrix, with endocervical and vaginal swabs being tested separately. Across all matrices, a total of 10 endogenous conditions and 6 exogenous substances were evaluated. Spiked samples were monitored for positivity and threshold inhibition. No interference was seen from the substance concentrations listed in Table 7.

	_	
1		Neisseria meningi Serogroup A
12		Neisseria meningi
4		Serogroup B
1		Neisseria meningi
3		Serogroup D
4	]/	Neisseria meningi Serogroup Y
2		Neisseria meningi
15		Serogroup W13
2		Neisseria cinere
2	\	Neisseria elonga
3	1 \	Neisseria flavesc
8		Neisseria lactam
		Neisseria muco
		Neisseria sicca
1		Neisseria subfla
51		Neisseria perfla
4		Neisseria polysacc
	12 4 1 3 4 2 15 2 2 3 3 8 1 2 3 8 1 1 1 1 51 4	12         4         1         3         4         2         15         2         15         2         3         4         2         15         2         3         8         1         51         4

**Table 6.** List of organisms used in the analytical specificity
 and cross-reactivity studies for the NeuMoDx CTNG Test. All organisms tested produced a negative result in absence of either CT or NG target. Additionally, the presence of these organisms at low levels did not affect the detection of CT or NG.

NeuMoDx CTNG Test - Endogenous/Exogenous						
Urine	3x LoD	10x LoD				
Blood, 7%	100%	100%				
Bilirubin, >1 mg/dL*	100%	100%				
Glucose, 1000 mg/dL	100%	100%				
Acidic Urine, pH 4	100%	100%				
Alkaline Urine, pH 9	100%	100%				
Protein (Albumin), 50 mg/mL	100%	100%				
PBMCs, 1x10 <sup>6</sup> cells/mL	100%	100%				
Talauma Davudan 19/	10.00/	100%				
Talcum Powder, 1%	100%	100%				
Swab	3x LoD	10x LoD				
Swab Blood, 10%	100% 3x LoD 100%	100% 10x LoD 100%				
Swab Blood, 10% Mucin, 12.9 mg/mL**	100% 3x LoD 100% 100%	100% 10x LoD 100% 100%				
Falcum Powder, 1%SwabBlood, 10%Mucin, 12.9 mg/mL**PBMCs, 1x105	100% 3x LoD 100% 100%	100% 10x LoD 100% 100%				
Faicum Powder, 1%SwabBlood, 10%Mucin, 12.9 mg/mL**PBMCs, 1x105Monistat 1, 21.2 mg/mL	100% 3x LoD 100% 100% 100%	100% 10x LoD 100% 100% 100%				
SwabBlood, 10%Mucin, 12.9 mg/mL**PBMCs, 1x105Monistat 1, 21.2 mg/mLKY Jelly Lubricant, 42.2 mg/mL	100% 3x LoD 100% 100% 100% 100%	100% 10x LoD 100% 100% 100% 100%				
Taicum Powder, 1%SwabBlood, 10%Mucin, 12.9 mg/mL**PBMCs, 1x105Monistat 1, 21.2 mg/mLKY Jelly Lubricant, 42.2 mg/mLYeast Gard Douche, 31.1 mg/mL	100%         3x LoD         100%         100%         100%         100%         100%         100%	100% 10x LoD 100% 100% 100% 100% 100%				
Taicum Powder, 1%SwabBlood, 10%Mucin, 12.9 mg/mL**PBMCs, 1x10 <sup>5</sup> Monistat 1, 21.2 mg/mLKY Jelly Lubricant, 42.2 mg/mLYeast Gard Douche, 31.1 mg/mLProgesterone, 5.9 mg/mL	100%         3x LoD         100%         100%         100%         100%         100%         100%         100%         100%         100%	100% 10x LoD 100% 100% 100% 100% 100% 100%				
Taicum Powder, 1%SwabBlood, 10%Mucin, 12.9 mg/mL**PBMCs, 1x10 <sup>5</sup> Monistat 1, 21.2 mg/mLKY Jelly Lubricant, 42.2 mg/mLYeast Gard Douche, 31.1 mg/mLProgesterone, 5.9 mg/mLVagisil, 7.1 mg/mL	3x LoD         100%         100%         100%         100%         100%         100%         100%         100%         100%	100% 10× LoD 100% 100% 100% 100% 100% 100%				

 
 Table 7.
 The NeuMoDx CTNG Test was performed with CTNG
 positive urine and swab matrix prepared by spiking negative matrix with CTNG targets at 3x LoD and 10x LoD in the presence of the listed interfering substances. Urine, n=3, Swab n=4.

\*Disease state clinical specimen \*\*Dosed from a 0.8% stock

### **Cross-Contamination**

The vulnerability of the NeuMoDx CTNG Test in conditions inclined to high levels of contamination was initially evaluated for cross contamination by performing a full run with alternating very high CTNG-positive urine and CTNG-negative urine in a "checkerboard" pattern. The possibility for carryover contamination was evaluated by running a full run of very high positive sample followed by a full run of negative samples. Both of these conditions were tested on two different NeuMoDx 288 Molecular Systems. For all valid results, no amplification was seen in any negative sample and all positive samples were correctly called positive for both targets.

### Specimen **Stability**

The NeuMoDx CTNG Test in Swab was evaluated for successful stability of specimens stored in normal refrigerator conditions with CTNG-negative clinical vaginal and endocervical swabs spiked with purified Chlamydia trachomatis (CT) elementary bodies (EB) and cultured *Neisseria gonorrhoeae* (NG) cells and stored between 2-8°C. At 30 days, 100% concordance with time O results was seen. Studies to assess specimen stability beyond 30 days are currently ongoing.

# **Turnaround Time**

The NeuMoDx CTNG Test on the NeuMoDx Molecular System was assessed for efficiency in turnaround time, which is defined as the time in minutes for samples to return a 'Completed' result after loading on the System for sample processing. Eight samples of CTNG-dosed urine were run on two different NeuMoDx Molecular Systems. Elapsed time was measured by checking the system reported start and end times, as well as a manual measurement with a handheld timer. All samples were completed in under an hour on both machines.

# NeuMoDx Molecular System Streamlined Testing



		1	NeuMoDx CTNG	Test in Urine-C	Contamination		
System	Level	<b>#POS Results</b>	<b>#NEG Results</b>	<b>#IND Results</b>	<b>#UNR Results</b>	Total #Valid Results	Total # Resul
			CROSS	<b>CONTAMINAT</b>	ION		
N000001	High Positive	24	0	0	0	24	24
NUUUUUI	Negative	0	23	0	0	23	23
N000004	High Positive	24	0	0	0	24	24
	Negative	0	23	0	1	23	24
TOTAL	High Positive	48	0	0	0	48	48
	Negative	0	46	0	1	46	47
			CARRYOV	ER CONTAMIN	ATION		
N00001	High Positive	46	0	1	0	46	47
NUUUUUI	Negative	0	47	0	0	47	47
	High Positive	48	0	0	0	48	48
N000004	Negative	0	46	2	0	46	48
TOTAL	High Positive	94	0	1	0	94	95
IUIAL	Negative	0	93	2	0	93	95

Table 8. Contamination data for both cross contamination and carryover contamination. All valid results show complete absence of contamination in either condition.

Specimen Stability of the NeuMoDx CTNG Test-Vaginal Swab							
		Confirmed CT	Positive Samples	Confirmed NG P	ositive Samples	Confirmed Neg	jative Samples
Timepoint		<b>#Positive</b>	#Negative	#Positive	#Negative	#Positive	#Negative
1	Time O	32/32	0/32	32/32	0/32	0/16	16/16
2	4 Days	32/32	0/32	32/32	0/32	0/15	15/15
3	7 Days	32/32	0/32	30/30	0/30	0/16	16/16
4	14 Days	32/32	0/32	32/32	0/32	0/16	16/16
5	30 Days	32/32	0/32	32/32	0/32	0/16	16/16
% Conco	ordance	10	0%	100	)%	100	0%
_	_	<u>Cuasius au C</u>			t Endersmiteel S		
		Specimen a	Desitive Semples	Confirmed NG P	st-Endocervical S	Confirmed Neg	nativo Samalor
Timepoint		#Positive	#Negative	#Positive	#Negative	#Positive	#Negative
1	Time O	32/32	0/32	32/32	0/32	0/16	16/16
2	4 Days	32/32	0/32	32/32	0/32	0/16	16/16
3	7 Days	32/32	0/32	32/32	0/32	0/16	16/16
4	14 Days	32/32	0/32	32/32	0/32	0/16	16/16
5	30 Days	32/32	0/32	32/32	0/32	0/16	16/16
% Conco	ordance	10	0%	100%		100%	
_	_	Snec	imen Stability of	the NeuMoDy CT	NG Test-I Irine	_	_
		Confirmed CT F	Positive Samples	Confirmed NG P	ositive Samples	Confirmed Neo	ative Samples
Timepoint		#Positive	#Negative	#Positive	#Negative	#Positive	#Negative
1	Time O	10/10	10/10	10/10	10/10	0/10	10/10
2	4 Days	10/10	10/10	10/10	10/10	0/10	10/10
3	7 Days	10/10	10/10	10/10	10/10	0/10	10/10
0		•	•	•			

e 9. The stability of the MoDx CTNG Test in Swab valuated at increasing points to assess test acy after sample storage °C. Invalid samples rerun until a valid result otained, and all results processed with the most -date assay definition neters.

**10.** The stability of the MoDx CTNG Test in Swab valuated at increasing points to assess test racy after sample storage 8°C. Invalid samples e rerun until a valid result obtained, and all results processed with the most -date assay definition neters.

**e 11.** The stability of the MoDx CTNG Test in Urine evaluated at increasing points to assess test uracy after sample storage

Turnaround Time of the NeuMoDx CTNG Test								
Machine	Method	Start Time	End Time	Elapsed Time	<b>Overall Results</b>	Replicates		
N00002	System Log	10:08 am	11:07 am	59 min.	Desitive	n=8		
N000002	Timer	0:00	0:58:40	~59 min.	Positive			
N000005	System Log	10:09 am	11:07 am	58 min.	Desitive			
11000005	Timer	0:00	0:58:42	~59 min.	FOSITIVE	n=8		

Table 12. Time elapsed was calculated from exported results on both machines. A manual measurement with a laboratory timer was also tracked. Runs on both machines completely processed 8 samples in less than 60 minutes.

### FEATURES

- Integrated Operation: Automates all steps of molecu diagnostics starting from raw clinical specimens to providing real-time PCR results in a fully automated
- **True Random Access:** Ability to mix specimen types and
- High Throughput: ~300 DNA tests in an 8 hour shift for the N288, ~150 DNA tests in an 8 hour shift for the N96
- Fast Time to First Results: ~60 min Continuous Loading: Specimens and Reagents can be
- loaded/unloaded at any time Large Walk-Away Window: Up to 288 samples for the
- N288. 96 samples for the N96 Seamless On Demand Operation: Automated inventory
- management of consumables and reagents Long In-Use shelf life: On-board room temperature
- Real-time PCR: Five-color fluorescence detection offers real-time PCR multiplexing ability





enables eas

cartridges capable of real-time PCR.

## Method Comparison Studies

The comparative performance of the NeuMoDx CTNG Test to an outside laboratory's FDA-cleared CTNG molecular diagnostic test was assessed for accuracy in both urine and swab matrices. The outside laboratory's diagnosti test was performed on individual patient specimens, and the residuals of these samples were then sent to NeuMoDx Molecular. Residual samples were analyzed by the NeuMoDx CTNG Test on two different NeuMoDx Molecular Systems. Results from NeuMoDx testing were compared to outside laboratory results to asses the NeuMoDx CTNG Test's sensitivity and specificity In urine, the overall sensitivity for the Chlamydia target was found to be 96.7% with

a specificity of 99.7%. The Gonorrhea target in urine was found to have a sensitivity of 98.1% and a specificity of 99.7%.

Swab samples were analyzed separately from urine, as well as being separated by type (endocervical vs. vaginal swabs). The overall sensitivity for both types of swab at the Chlamydia target was 100%, with a specificity of 99.6%. The Gonorrhea target saw a sensitivity of 100% and a specificity of 98.7% for all swab types combined. The accompanying tables show a

Method Comparison of the NeuMoDx CTNG Test-Urine						
Chlamydia		NeuMoDxTest				
		POS	NEG	TOTAL		
Reference	POS	87	3	90		
	NEG	1	297	298		
lest	TOTAL	88	300	388		
Sensitivity = 96.7% 95% CI (90.7% - 98.9%)						

Specificity = 99.7% 95% CI (98.1% - 99.9%)

Teet	NEG	1	335	336			
lest	TOTAL	52	336	338			
Sensitivity = 98.1% 95% CI (88.4% - 99.9%)							
Specificity = 99.7% 95% CI (98.1% - 100%)							

Table 13. Truth tables for the results of the NeuMoDx CTNG Test in comparison to FDA approved tests from reference labs as CT results and NG results.

Method Comparison of the NeuMoDx CTNG Test-Swab, Overall						
Chlannelia		NeuMoDxTest				
Chiam	iyala	POS	NEG	TOTAL		
Reference Test	POS	62*	0	62		
	NEG	1	263	264		
	TOTAL	63	263	326		
Sensitivity = 100% 95% CI (92.7-100%)						
Snoo	:f:_:t 00	6% 05% C	1/07 6 100	0/)		

Specificity = 99.6% 95% CI (97.6-100%) 24 of these specimens were a part of a contrived panel

targets (overall results for both swab types combined

Method Comparison of the NeuMoDx CTNG Test- Method Comparison of the NeuMoDx CTNG Test-

Endocervical Swabs						
Chlamydia		NeuMoDxTest				
		NEG	TOTAL			
POS	37*	0	37			
NEG	1	130	131			
TOTAL	38	130	168			
Sensitivity = 100% 95% CI (88.3-100%)						
Specificity = 99.2% 95% CI (95.2-100%)						
	Endo ydia POS NEG TOTAL itivity = 100	Endocervical SwhydiaNPOS37*POS37*NEG1TOTAL38itivity = 100% 95% CIificity = 99.2% 95% C	Endocervical Swabs           NeuMoDxTes           pydia         POS         NEG           POS         37*         0           NEG         1         130           TOTAL         38         130           itivity = 100%         95% CI (88.3-100)           ificity = 99.2%         95% CI (95.2-100)			

\*12 of these specimens were a part of a contrived panel.

argets (results for endocervical swabs only).					
Method Comparison of the NeuMoDx CTNG Test-Vaginal Swabs					
NeuMoDxTest				st	
Chlamydia		POS	NEG	TOTAL	
Reference Test	POS	25*	0	25	
	NEG	0	133	133	
	TOTAL	25	0	158	
Sensitivity = 100% 95% CI (83.4-100%)					
Specificity = 100% 95% CI (96.5-100%)					

	S	pecificity	/ =	98	.1%	95%	CI	(92)	.8-9
*48	of these	specimens	s we	ere a	part	of a cc	ntri∖	ved p	anel.

\*12 of these specimens were a part of a contrived panel. further breakdown of sensitivity Table 16. Truth tables for the results of the NeuMoDx CTNG Test in comparison to FDA approved tests from reference labs for CT and NG and specificity by swab type. targets (results for vaginal swabs only).

# **CONCLUSIONS & ACKNOWLEDGMENTS**

The NeuMoDx CTNG Test provides an extremely easy to use, reliable, low-cost, and rapid method for detection and differentiation of CT and NG directly from urine and swab specimens. For Research Use Only. Not for use in diagnostic procedures. | The authors gratefully acknowledge the help and support provided by all members of the NeuMoDx team. We would also like to thank our collaborators: Beaumont Health, Indiana University School of Medicine, Tampa General Hospital and Laboratory Alliance for their invaluable assistance in providing access to clinical specimens for testing and evaluation.

Urine	Method Co	mparison of	the NeuMo	Dx CTNG T	
	Comonstan		NeuMoDxTe		
OTAL	Gonor	rnea	POS	NEG	
90	Defense	POS	51	1	
298	Test	NEG	1	335	
388		TOTAL	52	336	
6 <b>)</b>	Sensiti	vity = 98.1%	% 95% CI (	88.4% - 99	
			NeuMoDxTes           POS         NEG           51         1           1         335           52         336           .1%         95% CI (88.4% - 99		

Method Comparison of the NeuMoDx CTNG Test-Swab, Overall					
Gonorrhea		NeuMoDxTest			
		POS	NEG	TOTAL	
Reference Test	POS	103*	0	103	
	NEG	3	220	220	
	TOTAL	106	220	326	
Sensitivity = 100% 95% CI (95.5-100%)					
Specificity = 98.7% 95% CI (95.8-99.7%)					

\*94 of these specimens were a part of a contrived panel.

**Table 14.** Truth tables for the results of the NeuMoDx CTNG Test in comparison to FDA approved tests from reference labs for CT and NG

Endocervical Swabs					
Gonorrhea		NeuMoDxTest			
		POS	NEG	TOTAL	
Reference Test	POS	52*	0	52	
	NEG	1	115	116	
	TOTAL	53	115	168	
Sensitivity = 100% 95% CI (91.4-10%)					
Specificity = 99.1% 95% CI (94.6-100%)					

\*46 of these specimens were a part of a contrived panel.

Table 15. Truth tables for the results of the NeuMoDx CTNG Test in comparison to FDA approved tests from reference labs for CT and NG

Method Comparison of the NeuMoDx CTNG Test-Vaginal Swabs					
Gonorrhea		NeuMoDxTest			
		POS	NEG	TOTAL	
Reference Test	POS	51*	0	51	
	NEG	2	105	107	
	TOTAL	53	105	158	
Sensitivity = 100% 95% CI (91.3-100%)					
Specificity = 98.1% 95% CI (92.8-99.7%)					